






Guide to Disinfectant Wipes

This reference guide is to assist when determining the best alternate cleaning wipe to use on a surface.

Type of Wipe	Approved Alternate cleaning wipe	Type of surface examples
 <p>Oxivir TB Clinical use (Small Tubs)</p> <p>Large tubs of Oxivir TB are strictly reserved for CSS use only.</p>	 <p>Accel Intervention Peroxigard</p>  <p>PREempt Accel Prevention</p>	<ul style="list-style-type: none"> ✓ Non porous surfaces ✓ Mobile equipment ✓ High touch surfaces ✓ Patient/non patient areas <p>Not recommended on</p> <ul style="list-style-type: none"> ✗ LCD screens/monitors ✗ Plexiglas ✗ Keyboards/computer screens
 <p>Cavi Wipes</p>	 <p>Micro-Kleen 3 Minuteman NEXGEN</p> <p>Pre COVID if your unit did not utilize the CAVI wipes continue with your established process</p>	<ul style="list-style-type: none"> ✓ LCD screens/monitors ✓ Specialty medical equipment ✓ Porous surfaces ✓ Plexiglas ✓ Plasma Screens/TV ✓ Lenses